

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA


Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision[®] USA



C400 MNRead Chart Black Text on White Background – In English

- The MNREAD Acuity Charts are continuous text reading acuity charts for normal and low vision.
- Often times traditional vision charts only measure general parts of normal vision, such as acuity and contrast.
- However, many times everyday life scenarios place these two aspects in unique situations.
- For example, acuity is very important in vision. Testing acuity normally involves a standard eye test. However, when reading the newspaper, the eyes are challenged by small print, continuous text, and several close contact lines.
- This basic living situation may cause problems for patients with scotomata and may not be detected on conventional acuity charts. Continuous text charts give one more possibility during vision testing.
- The research was funded and approved by the NIH.

Optotype	:	SLOAN
MN Read Charts	:	MN Read Chart 1-W, MN Read Chart 1-W Spanish, MN Read Chart 2-W, MN Read Chart 2-W Spanish, MN Read Chart 3-W, MN Read Chart 4-W, MN Read Chart 5-W
Chart Distance	:	40 cm / 16 in
Utilize With a Cabinet	:	No
Optotype Arrangement	:	Logarithmic Progression – Not ETDRS
Utilize with Which Cabinet	:	None
Patient Type	:	General / Literate
Contrast Level	:	High
Optotype	:	Text – Reading
ETDRS Format	:	No
Mass VAT Chart	:	No
Translucent Chart	:	No
Chart Size (cm)	:	30 x 46 cm
Cord	:	No
Chart Size (Inches)	:	12 x 18 in
Chart Notations – Snellen	:	20/400 to 20/6
Chart Notations – LogMAR	:	1.3 to – 0.5