

**Regd. Office :**

A- 6, NAND Complex,  
Near AIMS Oxygen, Juna Padra Road,  
Baroda 390 020 - INDIA


**Marketing / Customer Support :**

5-6, 3rd Floor, Mrudula Sadan,  
Pratap Road, Raopura,  
Baroda 390 001 - INDIA

**JUTRON Vision**

enhancing vision ~ enriching life

# Precision Vision<sup>®</sup> USA



## # 2126 Rabin Contrast Sensitivity Test For 914 Small Cabinet

- Contrast sensitivity test comparable to the Small Letter Contrast Test.
- Detects subtle effects of defocus, low luminance and glare.
- Log CS increases 0.25 steps per row (0.05 log CS/letter).
- Each line on the Rabin Contrast Sensitivity Test is a lower contrast than the one before it, using Log CS steps. This progression in contrasts helps to identify the patient's contrast sensitivity. Furthermore, the test helps to detect subtle effects such as defocus and low luminance and glare.