

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

112080 Vergence Facility Prism 12D/3D – Magnetic

These vergence prisms quickly screen for binocular vision. Use these to get the patient to cycle between the two prisms up to 15 cycles per minute. Use at 16" (40 cm) with a 20/30 near point line. Prisms are 12D base-out and 3D base-in and 1.5" high (38mm).

The Vergence Facility Prism is available in :

Magnetic which allows the prism to be used as either mounted or loose