

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



130001 3M PRESS ON PRISM

3M™ Press-On™ Prisms, also known as Fresnel Prisms, provide a system for creating prismatic corrections.

3M™ Press-On™ Prisms offer the eye care professional and the patient several benefits:

- A simple, therapeutic, inexpensive, way to correct several visual disorders.
- Provide an immediate correction.
- They are a more comfortable, more cosmetically appealing treatment for strabismus and diplopia than conventional prisms.
- They add NO noticeable weight or thickness to the spectacle.
- Press-on prism material is flexible static vinyl that can easily be cut to shape with scissors to determine the acceptance of a proposed corrective prescription.
- Adheres with just water to existing lenses, yet it can be easily repositioned.
- Apply a Press-On prism to the back surface of one or both lenses of the eyeglasses with just water.
- It can be applied to the entire lens or to any region of the lens.
- A full range of 17 powers from 1.00 to 40.00 prism diopters allows you to broaden therapeutic uses & provide additional treatment options in your practice.

One lens per box.