

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



GENERAL INSTRUCTIONS

1. If prescribed, take this test with reading glasses on, and view at 11 inches.
2. **Cover one eye** with your hand and look at the white solid square in the center.
3. Test the opposite eye in a similar manner.
4. If any lines or squares appear wavy, distorted or otherwise irregular or if a blind area appears during the testing of either eye, consult a licensed eye specialist.

COMPLIMENTS OF
Good-Lite
1155 Jansen Farm Drive • Elgin, IL 60123 USA
Tel 847-841-1145 • www.good-lite.com

181700 Amsler Giveaways

- Each sheet is intended to be given to a patient to be used at home. Another option would be to give the patient an Amsler Viewer which has the added advantage that it can only be used monocularly.
- The grid is a black background with white lines. Sold as a package of 250 sheets.