

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA


Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision[®] USA



2110 Original Series Sloan Letter ETDRS (Chart R)

Designed for use in clinical studies and low vision evaluations where accurate, repeatable visual acuity measurements are required.

Can be wall-mounted.

ETDRS charts use a geometric progression from line to line, making only the optotype size a variable. Lines and optotypes have been calibrated to be equally difficult. ETDRS chart is reliable for accurate acuity testing in clinics, research, and studies world wide.