

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision[®] USA

Row	<small>RABIN CONTRAST SENSITIVITY TEST DEVELOPED BY JEFF RABIN, O.D., Ph.D. ETC IS GLASSMAY (PNEI) FOR TESTING AT 4 METERS (13 FEET) LETTER SIZE 20M, 4.0D, 10/16 (8/16, 20/70) SQUARE</small>					log CS
1	Z	R	K	D	C	0.25
2	D	N	C	H	V	0.50
3	C	D	H	N	R	0.75
4	R	V	Z	O	S	1.00
5	O	S	D	V	Z	1.25
6	N	O	Z	C	D	1.50
7	R	D	N	S	K	1.75
8	O	K	S	V	Z	2.00

2128-1 Rabin Contrast Sensitivity Test For #2425 Large Cabinet

- Detects subtle effects of defocus, low luminance and glare.
- Log CS increases 0.25 steps per row (0.05 Log CS per Letter).
- Must be used with Cat. No. 2425 Large ETDRS Illuminator Cabinet.
- Rabin Contrast Sensitivity Test uses several different levels of contrast to find the patient's true sensitivity. The line progression helps to detect early signs such as defocus, glare and low luminance problems. Because of the precision of each contrast level, the large Illuminator Cabinet (2425) must be used. Without the use of the cabinet, the contrast levels will not be accurate.