

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA


Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision[®] USA



2424 Large mesopic Filter for 2425 ETDRS Illuminator Cabinet

Thin and light design

Testing for multiple vision problems while performing routine exams is practical. Acuity, Low Contrast and Low Vision testing are important parts of normal vision and should be checked and recorded often. Placed in front of a vision test chart, this filter reduces normal cabinet light level down to 3 cd/m² to challenge the patient's vision at low light conditions. It is an easy and quick way to test low vision and low lighting problems and does not require moving to different rooms, calculating, or reorganizing the office. Mesopic filter is light, easy to use, and creates an accurate low vision test.