

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

346707 Handheld Bagolini Striated Lens Set With Penlights

- These finely striated clear lenses are used for evaluating retinal correspondence in patients with eye deviations and in the analysis of anomalous correspondence, suppression, etc. Consists of 2 circular lenses mounted in a flipper.
- When you use the penlight, a patient with no suppression will see two diagonal lines crossing above, below, or at the light source. This can be accomplished at a near or far distance. If suppression is present, all or part of one line will not be seen. If the diagonal lines cross at the source when the cover test indicated an ocular deviation, the patient has harmonious, abnormal retinal correspondence.
- Penlight has a pupillometer ranging from 2 to 9.