Regd. Office:
A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support: 5-6, 3rd Floor, Mrudula Sadan, Pratap Road, Raopura, Baroda 390 001 - INDIA







## # 494322 Color Loss Simulating Flipper

This flipper uses frequency isolation to simulate Color Vision, loss. This product is used to help patients understand what color deficiency is, and how it affects vision.

This flipper uses frequency isolation to simulate Color Vision, loss. This product is used to help patients understand what color deficiency is, and how it affects vision.