

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



607700 Optokinetic Flag – Black and White 1.5" Stripes

This optokinetic flag is used to produce optokinetic nystagmus (OKN) physiologic jerk in a patient's eyes. This is elicited when the patient attempts to fixate on moving, alternative, Contrasting squares that are moved in the field of vision.

This Optokinetic flag is available with 1½" (4 cm) segments. Comes with a guide and is 27" (68 cm) long.