

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



484050-1 Jackson Cross Cylinder 0.25

These Jackson Cross Cylinders are used to determine the orientation and power of the astigmatism correction and the axis. The lenses are 35mm with axis markings and a handle that is 70mm long.