

**Regd. Office :**

A- 6, NAND Complex,  
Near AIMS Oxygen, Juna Padra Road,  
Baroda 390 020 - INDIA

**Marketing / Customer Support :**

5-6, 3rd Floor, Mrudula Sadan,  
Pratap Road, Raopura,  
Baroda 390 001 - INDIA

**JUTRON Vision**

enhancing vision ~ enriching life



## # BSDSH1 Dual Horizontal Slide Holder

The dual horizontal slide holder is an economical option that allows the patient to train in free space, shift accommodation through the paddles and back to the target. By turning the paddles around you can easily shift from base-in to base-out demand, as well as change the 3D target from coming towards you to going away. Perform jump duction training with two different base-in demands or two base-out demands or have top base-in and bottom base-out and "jump" back and forth between the two images.

The clear fixation paddle is 12" x 14" made of clear polycarbonate. It utilizes (2) horizontal slide holders that are attached to the paddle via Velcro® connections. These slide holders fit Stereo Optical Vectograms® as well as Bernell Tranaglyphs™, as well as Bernell Clings.

**Includes:** (1) Clear Acrylic Fixation Paddle with Handle & (2) Clear Horizontal Slide Holders with attachments

\*Does not include pictured Tranaglyph and Vectogram

Delhi : Chandigarh : Lucknow : Jaipur : Indore : Ahmedabad : Mumbai : Pune : Nagpur : Hyderabad : Bangalore : Chennai

Tel. : +91-99987-02020 • mail@jutronvision.com • www.JutronLowVision.com • www.JutronVision.com