


Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

C+ 3M Press-On Aspheric Lenses

3M™ Press-On™ Optics, also known as Fresnel Prisms, and Press-On Aspherics, provide a system for creating prismatic and spherical corrections. 3M™ Press-On™ Optics offer the eye care professional and the patient several benefits:

- ☐ A simple, therapeutic, inexpensive, way to correct several visual disorders.
- ☐ Provide an immediate correction.
- ☐ They are more comfortable, more cosmetically appealing treatment for strabismus and diplopia than conventional prisms.
- ☐ They add NO noticeable weight or thickness to the spectacle.
- ☐ Adheres with just water to existing lenses, yet it can be easily re-positioned.
- ☐ Apply a Press-On sphere or prism to the back surface of one or both lenses of the patient's eyeglasses with just water.
- ☐ A range of powers from -13.00 to +16.00 prism diopters allows you to broaden therapeutic uses & provide additional treatment options in your practice. Made of polyvinyl chloride