
Regd. Office:
A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support: 5-6, 3rd Floor, Mrudula Sadan, Pratap Road, Raopura, Baroda 390 001 - INDIA







GNTCTC2 Gottlieb Near Training Cards: Column of Circles™

This advanced VT procedure is used to fine tune accommodation and convergence/divergence. The circles appear at a great deal of variation in depth. Letters in the circles given even more variations in depth. The patient's attention is held for a long time to fine tune their accommodative and convergen/divergence systems.

Includes:

- 1 Laminated Cardstock Card
- 1 Acetate Card
- Instruction Manual

Size: 8.5" x 11"