

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



JUPOK Jazz Ultra Pocket Ophthalmoscope

Stylish Customisable Pocket Ophthalmoscope for Ward and GP use

The Jazz Ultra has a more compact design than the Jazz, with a smaller and more portable handle.

Full diagnostic capacity in a compact instrument make the Pocket Ophthalmoscope ideal for the Practitioner on the move. An ideal choice for domiciliary and ward use or as a backup instrument.

Our stylish Jazz Ophthalmoscope will fit in any environment. You can change the colour of your set by adding coloured rings that can suit any mood.

Four beams for increased control over diagnostic functionality.

Additional blue filter that allows assessment of corneal scratches.

2.8V model accepts AA Batteries. Jazz Ultra Handle not compatible with other ophthalmoscopes.

Optical Excellence, Manufactured in the UK

Contents:

Jazz Ultra LED Pocket ophthalmoscope

Jazz Ultra Handle

Protective case

Interchangeable Colour Rings

Features:

Our most compact pocket ophthalmoscope.

Delhi : Chandigarh : Lucknow : Jaipur : Indore : Ahmedabad : Mumbai : Pune : Nagpur : Hyderabad : Bangalore : Chennai

Tel. : +91-99987-02020 • mail@jutronvision.com • www.JutronLowVision.com • www.JutronVision.com

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

Slimline design fits neatly in your pocket.

Lens powers range from +20D to -20D.

Interchangeable Colour Rings allows Customisation.

Four diagnostic beams for full diagnostic capability.

Wide Angle – for general examination.

Macula – for more detailed views.

Semi Circle – combination of depth perception and field of view

Fixation Cross – projects a graticule onto the retina for assessment of eccentric fixation

Green, Red Free, Filter.

Blue Filter – to be used in combination with Fluorescein for assessing corneal scratches.