

**Regd. Office :**

A- 6, NAND Complex,  
Near AIMS Oxygen, Juna Padra Road,  
Baroda 390 020 - INDIA


**Marketing / Customer Support :**

5-6, 3rd Floor, Mrudula Sadan,  
Pratap Road, Raopura,  
Baroda 390 001 - INDIA

**JUTRON Vision**

enhancing vision ~ enriching life

# Precision Vision<sup>®</sup> USA



## # 2004 HOTV Intermediate Vision Chart 40cm

Near Vision Testing at 40cm: 8M (20/400) to .20M (20/10)

Memorization can alter testing results and create an inaccurate record for the patient. This near vision chart features multiple groups for repetitive testing help minimize memorization. This helps to create a more accurate and reliable vision chart that yields repeatable results.

HOTV: Standard Near 40cm, MassVat