

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA


Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision[®] USA



2106-66 SLOAN ETDRS Intermediate Vision Chart (66cm)

- Intermediate at a Choice of 66cm Testing Distances.
- Chart with Cord for easily achieving proper testing distance.
- Intermediate Vision Testing at 66cm: 8M (20/250) to .20M (20/6.3)
- Memorization can alter testing results and create an inaccurate record for the patient. This near vision chart features multiple groups for repetitive testing help minimize memorization. This helps to create a more accurate and reliable vision chart that yields repeatable results.
- Each side is different to minimize memorization.
- Sloan ETDRS: Intermediate 66cm