

**Regd. Office :**

A- 6, NAND Complex,  
Near AIMS Oxygen, Juna Padra Road,  
Baroda 390 020 - INDIA


**Marketing / Customer Support :**

5-6, 3rd Floor, Mrudula Sadan,  
Pratap Road, Raopura,  
Baroda 390 001 - INDIA

**JUTRON Vision**

enhancing vision ~ enriching life

# Precision Vision<sup>®</sup> USA



## # 2211 Landolt C Intermediate Vision Chart 40cm

Memorization can alter testing results and create an inaccurate record for the patient. This near vision chart features multiple groups for repetitive testing help minimize memorization. This helps to create a more accurate and reliable vision chart that yields repeatable results.

Landolt C: Standard Near 40cm