

Regd. Office :


A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life



SOV2 Variable Quoits Vectogram

Variable Quoits Vectogram

- 8-1/4" x 5-1/2"
- Polarized vectograms ensure testing and training at actual distances. vectograms stabilize fusion and stereopsis, eliminate suppression and develop simultaneous vision.
- Used in monocular, binocular, pursuit and saccadic training, vectograms help correct and prevent anomalies involving projection and hand-eye coordination.

Ideal for use with Bernell Dual Polachrome™ Trainer

Polarized Goggle Included

Cleaning Instructions:

- **DO NOT IMMERSE THE VECTOGRAM PANELS IN WATER**
- To clean the vectogram panels and guide, use a soft, damp, lint-free cloth. Dampen cloth using glass cleaner or mild detergent/water.
- If panels are removed from guide during cleaning, replace the panels in the guide, placing the panel with the solid bar on top of the panel with the numbers/letters and ensure that the plastic portion of the guide is behind the panels

Storage Instructions:

- Store Vectogram in a dry place, away from direct sunlight.

Warranty: 1 year limited manufacturer warranty from the date of purchase