

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

251300 Lea Grating Acuity Test

This test at full contrast requires discrimination of the direction of long lines, which is a more demanding task than the detection acuity measured with Lea GRATINGS.

The test gratings are presented to children and adults at different distances to determine which distance the direction of lines can be perceived.

Knowing the distance and which grating was used, the grating acuity values can be read from the nomograms. Includes, nomograms, instructions and storage case.

4 disks 8" (20 cm) in diameter.