

**Regd. Office :**

A- 6, NAND Complex,  
Near AIMS Oxygen, Juna Padra Road,  
Baroda 390 020 - INDIA

**Marketing / Customer Support :**

5-6, 3rd Floor, Mrudula Sadan,  
Pratap Road, Raopura,  
Baroda 390 001 - INDIA


**JUTRON Vision**
*enhancing vision ~ enriching life*


Chart 1

Chart 2

## # 500330 Evans Letter Contrast Test (ELCT) – Set

The Evans Letter Contrast Test (ELCT) chart is the first standardized test for letter contrast. The ELCT chart is designed to be used with either the ESV3000 or ESC2000 providing consistent lighting for improved testing results. The ELCT chart presents Sloan Letter, optotypes line size of 20/630 with a testing distance of 1m. The varying contrast in sets of 3 letters of equal contrast level. All triplets are balanced for letter difficulty and differ in contrast by 0.15 log contrast step.

**ELCT Chart 1** contains 12 contrast levels ranging from standard 100% acuity to 1.65 log contrast or 2.24% contrast level. Studies show that 1.65 log contrast is the average letter contrast for the normal older population. (Weber Contrast calculation).

**ELCT Chart 2** presents 4 contrast levels from 1.8 to 2.25 log contrast.

Set includes both Chart 1 and Chart 2.

Dimensions: Height 24.5" x Width 25.5" (62 cm x 65 cm).