

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

657500 WOLFF WAND FIXATION TARGETS WITH CASE

These Wolff Wands provide specialized fixation targets used in various procedures such as retinoscopy, cover tests, and ocular motility. Various techniques employ the wands to evaluate accommodative lag, anisometropia, near point convergence, tracking, peripheral awareness, and saccadic eye movements.

This basic Wolff Wand set contains two 14" (35.6 cm) steel rods: one with a 1/2" (13 mm) polished brass ball and one with a silver nickel-plated ball. The rods are made of aluminum for durability. The balls are polished to permit fixation on either the ball or the patient's image in the shiny ball.