

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life***# PocketOK****Pocket Ophthalmoscope**

Portable & Compact Pocket Ophthalmoscope for Ward and GP use

Full diagnostic capacity in a compact instrument make the Pocket Ophthalmoscope ideal for the Practitioner on the move. An ideal choice for domiciliary and ward use or as a backup instrument.

Five beams for increased control over diagnostic functionality.

2.8V model accepts AA Batteries. Pocket Handle not compatible with other Ophthalmoscopes.

Optical Excellence, Manufactured in the UK

Contents:

Pocket ophthalmoscope

Pocket Handle

Spare bulb

Instructions CD

Protective case

Features:

Slimline design fits neatly in your pocket.

Lens powers range from +20D to -20D.

Strong metal clip that can be attached to a pocket.

Five diagnostic beams for full diagnostic capability.

Wide Angle – for general examination.

Macula – for more detailed views.

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA

Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision*enhancing vision ~ enriching life*

Intermediate – providing easier access for non-dilated pupils and in peripheral examination. Ideal for paediatrics.

Slit – to determine retinal elevations and depressions. Used to access the anterior chamber depth.

Fixation Cross when used in combination with the Red Free Filter – projects a graticule onto the retina for assessment of eccentric fixation.

Swing over Green, Red Free, Filter.

SKU: PocketOK Categories: Examinati