

Regd. Office :

A- 6, NAND Complex,
Near AIMS Oxygen, Juna Padra Road,
Baroda 390 020 - INDIA


Marketing / Customer Support :

5-6, 3rd Floor, Mrudula Sadan,
Pratap Road, Raopura,
Baroda 390 001 - INDIA

JUTRON Vision

enhancing vision ~ enriching life

Precision Vision®



2103 Logarithmic SLOAN Visual Acuity Test

Two smaller groups 6M to 1.5M for repetitive testing (20/40, 6/12 to 20/10, 6/3 equiv.).

Can be wall-mounted.

Includes LogMar and Letter sizes on chart.

Logarithmic chart uses traditional Sloan optotype

These calibrated optotype

s blur at threshold and adhere to the Snellen 5×5 grid. Sloan chart uses two smaller groups for repetitive testing and minimizes memorization.

Logarithmic Sloan acuity test is a great screener for literate adults and children.